

CLEAN DOSE PRESOAK TABLET

Enzymatic formula breaks down
protein and starch films



Multi-use and effective

- Presoaking flatware results in consistently clean and spotless results with every wash cycle
- Saves time and labor by helping to reduce need for rewash

User-friendly design

- Safer product handling
- Requires no PPE

Easy on the environment

- Concentrated solid tablets contribute to reduced water, energy, and plastic waste

The tablet's solid chemistry
is formulated to save on
packaging.

99.9%
reduction in plastic waste¹

1 Reduction in plastic waste is an approximation assuming use of one case of product per month over 12 months, compared with the equivalent use of Ecolab liquid products, which is enabled by the smaller, lighter package size for solid tablets.

DIRECTIONS FOR USE

1. Tear open a packet and drop the tablet into a presoak tub. Wash hands after handling product.
2. Fill up the tub with warm water (110° – 120° F) and allow the tablet to dissolve.
3. Place dirty flatware into the solution and allow to soak.

CLEAN DOSE PRESOAK TABLET	
PACK SIZE	48 ct
DILUTION RATIO	1 tablet/tub
SWSH #	6102782
DIST #	



We thought small so you can think big.

**Ready to learn more about our full tablet cleaning program?
Visit swsh.com/cleandose or scan the QR code.**

The experienced and dedicated Swisher sales team provides you personalized support to help you run a clean and safe operation.

Call **1 800 360 SWSH** or visit **[SWSH.com](https://swsh.com)** for more information.

1 Ecolab Place St. Paul, MN 55102

©2025 Ecolab USA Inc. All rights reserved. 62074/0400/0225

